

# Simple Harmonic Motion

Lots of things vibrate or oscillate. A vibrating tuning fork, a moving child's playground swing, and the speaker in a headphone are all examples of physical vibrations. There are also electrical and acoustical vibrations, such as radio signals and the sound you get when blowing across the top of an open bottle.

One simple system that vibrates is a mass hanging from a spring. The force applied by an ideal spring is proportional to how much it is stretched or compressed. Given this force behavior, the up and down motion of the mass is called *simple harmonic* and the position can be modeled with

$$y = A\sin(2\pi ft + \phi)$$

In this equation,  $y$  is the vertical displacement from the equilibrium position,  $A$  is the amplitude of the motion,  $f$  is the frequency of the oscillation,  $t$  is the time, and  $\phi$  is a phase constant. This experiment will clarify each of these terms.




Figure 1

## OBJECTIVES

- Measure the position and velocity as a function of time for an oscillating mass and spring system.
- Determine the amplitude, period, and phase constant of the observed simple harmonic motion.
- Compare the observed motion of a mass and spring system to a mathematical model of simple harmonic motion.

## MATERIALS


computer  
Vernier computer interface  
Logger Pro  
Vernier Motion Detector  
200 g and 300 g masses

ring stand, rod, and right-angle clamp  
spring, with a spring constant of  
approximately 15 N/m  
twist ties  
wire basket

## PRELIMINARY QUESTIONS

1. Attach the 200 g mass to the spring and hold the free end of the spring in your hand, so the mass and spring hang down with the mass at rest. Lift the mass about 10 cm and release. Observe the motion. Sketch a graph of position vs. time for the mass.
2. Just below the graph of position vs. time, and using the same length time scale, sketch a graph of velocity vs. time for the mass.

## PROCEDURE

1. Attach the spring to a horizontal rod connected to the ring stand and hang the mass from the spring, as shown in Figure 1. Securely fasten the 200 g mass to the spring and the spring to the rod, using twist ties so the mass cannot fall. Adjust the height of the mass so that the bottom of the mass is about 33–55 cm from the table top or floor.
2. Set the Motion Detector sensitivity switch to Ball/Walk. Connect the Motion Detector to a digital (DIG) port of the interface. 
3. Place the Motion Detector below the mass. No objects should be near the path between the detector and mass, such as a table edge. Place the wire basket over the Motion Detector to protect it.
4. Open the file “15 Simple Harmonic Motion” from the *Physics with Vernier* folder.
5. Make a preliminary run to verify things are set up correctly. Lift the mass upward a few centimeters and release. The mass should oscillate along a vertical line only. Click  to begin data collection.
6. When data collection is complete, the position graph should show a clean sinusoidal curve. If it has flat regions or spikes, reposition the Motion Detector and try again.
7. Compare the position graph to your sketched prediction in the Preliminary Questions. How are the graphs similar? How are they different? Also, compare the velocity graph to your prediction.
8. Estimate the equilibrium position of the 200 g mass. Do this by allowing the mass to hang free and at rest. Click  to begin data collection. After collection stops, click Statistics,  $\frac{\sqrt{2}}{2}$ , to determine the average distance from the detector. Record this value as position ( $y_0$ ) for Run 1 in your data table.
9. Now lift the mass upward about 5 cm and release it. The mass should oscillate along a vertical line only. Click  to collect data. Examine the graphs. The pattern you are observing is characteristic of simple harmonic motion.
10. Use the position graph to measure the time interval between maximum positions. This is the period,  $T$ , of the motion. The frequency,  $f$ , is the reciprocal of the period,  $f = 1/T$ . Based on your period measurement, calculate the frequency. Record the period and frequency of this motion in your data table.

11. The amplitude,  $A$ , of simple harmonic motion is the maximum distance from the equilibrium position. Estimate values for the amplitude from your position graph. Enter the values in your data table. If you drag the mouse from a peak to an adjacent trough, *Logger Pro* will report the change in position over that region.
12. Repeat Steps 9–11 with the same 200 g mass, but with a larger amplitude than in the first run.
13. Change the mass to 300 g and repeat Steps 8–11. Use an amplitude of about 5 cm. Keep a good run made with this 300 g mass on the screen. You will use it for several of the Analysis questions.

### DATA TABLE

Run	Mass (g)	$y_0$ (m)	$A$ (m)	$T$ (s)	$f$ (Hz)
1					
2					
3					

Time (s)	Position (m)
when $v = 0$ _____	
when $v$ is maximum _____	

Fitted equation with parameters	
---------------------------------	--

### ANALYSIS

1. View the graphs of the last run. Compare the position *vs.* time and the velocity *vs.* time graphs. How are they the same? How are they different?
2. Click Examine, to use the Examine tool. Move the mouse cursor back and forth across the graph to view the data values for the last run on the screen. In your data table, record time and position values for when  $v = 0$ . Also record time and position values for a point when the velocity is greatest. Relative to the equilibrium position, where is the mass when the velocity is zero? Where is the mass when the velocity is greatest?
3. Does the frequency,  $f$ , appear to depend on the amplitude of the motion? Do you have enough data to draw a firm conclusion?
4. Does the frequency,  $f$ , appear to depend on the mass used? Did it change much in your tests?

5. You can compare your experimental data to the sinusoidal function model using the Model feature of Logger Pro. Try it with your 300 g data. The model equation in the introduction, which is similar to the one in many textbooks, gives the displacement from equilibrium. However, your Motion Detector reports the distance from the detector. To compare the model to your data, add the equilibrium distance to the model; that is, use

$$y = A \sin(2\pi ft + \phi) + y_0$$

where  $y_0$  represents the equilibrium distance. The phase parameter,  $\phi$ , is called the *phase constant* and is used to adjust the  $y$  value reported by the model at  $t = 0$  so that it matches your data.

- Click once on the position graph to select it.
  - Choose Model from the Analyze menu and select Latest.
  - Select the Sine function from the General Equation list.
  - The Sine equation is of the form  $y=A*\sin(Bt +C) + D$ . Compare this to the form of the equation above to match variables; e.g.,  $\phi$  corresponds to  $C$ , and  $2\pi f$  corresponds to  $B$ .
  - Adjust the values for  $A$ ,  $B$  and  $D$  to reflect your values for  $A$ ,  $f$  and  $y_0$ . You can either enter the values directly in the dialog box or you can use the up and down arrows to adjust the values.
  - The optimum value for  $\phi$  will be between  $0$  and  $2\pi$ . Find a value for  $\phi$  that makes the model come as close as possible to the data of your 300 g experiment. You may also want to adjust  $y_0$ ,  $A$ , and  $f$  to improve the fit. Write down the equation that best matches your data.
6. Does the model fit the data well? How can you tell?
7. Predict what would happen to the plot of the model if you doubled the parameter for  $A$  by sketching both the current model and the new model with doubled  $A$ . Now double the parameter for  $A$  in the model dialog box to compare to your prediction.
8. Similarly, predict how the model plot would change if you doubled  $f$ , and then check by modifying the model definition.

## EXTENSIONS

- Investigate how changing the spring amplitude changes the period of the motion. Take care not to use too large of an amplitude so that the mass does not come closer than 15 cm to the detector or fall from the spring.
- How will *damping* change the data? Tape an index card to the bottom of the mass and collect additional data. You may want to take data for more than 10 seconds. Does the model still fit well in this case?
- Do additional experiments to discover the relationship between the mass and the period of this motion.

# Vernier Lab Safety Instructions Disclaimer

**THIS IS AN EVALUATION COPY OF THE VERNIER STUDENT LAB.**

**This copy does not include:**

- **Safety information**
- **Essential instructor background information**
- **Directions for preparing solutions**
- **Important tips for successfully doing these labs**

The complete *Physics with Vernier* lab manual includes 35 labs and essential teacher information. The full lab book is available for purchase at: [www.vernier.com/pwv](http://www.vernier.com/pwv)



Vernier Software & Technology  
13979 S.W. Millikan Way • Beaverton, OR 97005-2886  
Toll Free (888) 837-6437 • (503) 277-2299 • FAX (503) 277-2440  
[info@vernier.com](mailto:info@vernier.com) • [www.vernier.com](http://www.vernier.com)